



1565 South Shields Drive, Waukegan, IL 60085
tel (847) 688-9968 fax (847) 688-9768
www.adimedical.com email: adimedical@msn.com

Product Specification

Synthetic (Latex Free) Tourniquets

Item Number 10910 / 10920 / 10930

A special mixture of polyisoprene rubber and other fine materials create a unique latex-free surface.

Features:

- ✓ Latex Free
- ✓ Characteristics, both physical and performance, similar to natural latex
- ✓ Easy to tie and untie
- ✓ Tear and slip resistant
- ✓ Suitable for medical kits
- ✓ Provide either flat or already rolled and banded
- ✓ Availability: ¾" and 1" widths
- ✓ Non-Sterile

Physical Properties *	Typical Results
Tensile Strength (mega pascal)	17Mpa Min
Elongation at break	800% Min
Deformation %	15%
Rigidity (A)	45 degrees

**Physical Properties are typical and may vary.*