



1565 South Shields Drive, Waukegan, IL 60085
tel (847) 688-9968 fax (847) 688-9768
www.adimedical.com email: adimedical@msn.com

Product Specification

Non-Woven Gauze Sponges

- A. Description:
- # 20024 Gauze sponge non-woven 2x2-4 ply, 5m/cs
 - # 20026 Gauze sponge non-woven 2x2-6 ply, 8m/cs
 - # 20034 Gauze sponge non-woven 3x3-4 ply, 4m/cs
 - # 20044 Gauze sponge non-woven 4x4-4 ply, 4m/cs
 - # 20826 Gauze sponge non-woven Slit 2x2-6 ply, 5m/cs
 - # 20847 Gauze sponge non-woven Slit 4x4-6 ply, 5m/cs
- B. Physical Properties:
- Rayon and viscos spun-laced non-woven fabric
 - White color
 - Non-Sterile
 - Durable and absorbent
- C. 510K number: exempt
- D. Latex free
- E. Made in China
- F. Medical device number: A877763